POLARIS

FR – NOTICE D'UTILISATION
EN - INSTRUCTIONS FOR USE
DE - GEBRAUCHSANWEISUNG
ES - INSTRUCCIONES DE USO
IT - ISTRUZIONI PER L'USO
PT - INSTRUÇÕES DE UTILIZAÇÃO
DA - BRUGSANVISNING
FI - KÄYTTÖOHJEET
NL - GEBRUIKSAANWIJZING
SV - BRUKSANVISNING
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
BG - ИНСТРУКЦИИ ЗА УПОТРЕБА
RO - INSTRUCȚIUNI DE UTILIZARE
ET - KASUTUSJUHEND
LV - LIETOŠANAS INSTRUKCIJA
HU - HASZNÁLATI UTASÍTÁS
SK - NÁVOD NA POUŽITIE
SL - NAVODILA ZA UPORABO
CS - NÁVOD K POUŽITÍ
MT - ISTRUZZJONIJIET GĦALL-UŻU
HR - UPUTE ZA KORIŠTENJE
GA - TREORACHA ÚSÁIDE
LT - NAUDOJIMO INSTRUKCIJOS
PL - INSTRUKCJA UŻYTKOWANIA

DESCRIPTION

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DESCRIPTION

POLARIS ultrasound or electrostimulation gel is a non-sterile, non-measuring Class I product according to Rule 1 of Annex VIII of the EU Medical Devices Regulation 2017/745 and is intended to facilitate the transmission of ultrasound between the ultrasound or electrostimulation machine and the patient. POLARIS gel is colorless, transparent, fragrance-free and paraben-free.

GUIDANCE FOR USE

Apply POLARIS gel on healthy skin and/or on the ultrasound probe in sufficient quantity for the examination to be performed. After each use, wipe off any gel remaining on the skin with a paper towel. Wipe the ultrasound probe with a damp cloth, then with a disinfectant approved by the probe manufacturer, between each patient.

STORAGE CONDITIONS

Keep away from direct sunlight, in a cool and dry place and at room temperature. POLARIS gel has a shelf life of 3 years. To be used in accordance with the French Directive No. 960479: can be kept for 24 hours after opening.

PRECAUTIONS

For external use only. Do not swallow. For medical use only. Any serious incident related to the device should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

DISPOSAL

Polyethylene bottle and cap, to be disposed of in accordance with current regulations.

NF EN ISO 15223-1 2021

	Titre	Description
	Medical device	Indicates that the item is a medical device
	CE marking	Indicates that the article complies with the requirements set by the regulations in force
	Batch code	Indicates the manufacturer's batch code in order to formally identify the batch
	Expiration date	Indicates the date after which the medical device should no longer be used
	Manufacturer	Indicates the manufacturer of the medical device
	Authorised representative in Switzerland	Indicates the authorised representative in Switzerland
	Consult the electronic instructions for use	Indicates the need for the user to consult the instructions for use

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