# SONECHO

FR – NOTICE D'UTILISATION  
EN - INSTRUCTIONS FOR USE  
DE - GEBRAUCHSANWEISUNG  
ES - INSTRUCCIONES DE USO  
IT - ISTRUZIONI PER L'USO  
PT - INSTRUÇÕES DE UTILIZAÇÃO  
DA - BRUGSANVISNING  
FI - KÄYTTÖOHJEET  
NL - GEBRUIKSAANWIJZING  
SV - BRUKSANVISNING  
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
BG - ИНСТРУКЦИИ ЗА УПОТРЕБА  
RO - INSTRUCȚIUNI DE UTILIZARE  
ET - KASUTUSJUHEND  
LV - LIETOŠANAS INSTRUKCIJA  
HU - HASZNÁLATI UTASÍTÁS  
SK - NÁVOD NA POUŽITIE  
SL - NAVODILA ZA UPORABO  
CS - NÁVOD K POUŽITÍ  
MT - ISTRUZZJONIJIET GĦALL-UŻU  
HR - UPUTE ZA KORIŠTENJE  
GA - TREORACHA ÚSÁIDE  
LT - NAUDOJIMO INSTRUKCIJOS  
PL - INSTRUKCJA UŻYTKOWANIA

## DESCRIPTION

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## DESCRIPTION

SONECHO ultrasound or electrostimulation gel is a non-sterile, non-measuring Class I product according to Rule 1 of Annex VIII of the EU Medical Devices Regulation 2017/745 and is intended to facilitate the transmission of ultrasound between the ultrasound or electrostimulation machine and the patient. SONECHO gel is colorless, transparent, fragrance-free and paraben-free.

## GUIDANCE FOR USE

Apply SONECHO gel on healthy skin and/or on the ultrasound probe in sufficient quantity for the examination to be performed. After each use, wipe off any gel remaining on the skin with a paper towel. Wipe the ultrasound probe with a damp cloth, then with a disinfectant approved by the probe manufacturer, between each patient.

## STORAGE CONDITIONS

Keep away from direct sunlight, in a cool and dry place and at room temperature. SONECHO gel has a shelf life of 3 years. To be used in accordance with the French Directive No. 960479: can be kept for 24 hours after opening.

## PRECAUTIONS

For external use only. Do not swallow. For medical use only. Any serious incident related to the device should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

## DISPOSAL

Polyethylene bottle and cap, to be disposed of in accordance with current regulations.

## NF EN ISO 15223-1 2021

|   | Titre                                       | Description   |
|---|---|---|
|  | Medical device                              | Indicates that the item is a medical device   |
|  | CE marking                                  | Indicates that the article complies with the requirements set by the regulations in force |
|  | Batch code                                  | Indicates the manufacturer's batch code in order to formally identify the batch           |
|  | Expiration date                             | Indicates the date after which the medical device should no longer be used                |
|  | Manufacturer                                | Indicates the manufacturer of the medical device  |
|  | Authorised representative in Switzerland    | Indicates the authorised representative in Switzerland                                    |
|  | Consult the electronic instructions for use | Indicates the need for the user to consult the instructions for use                       |

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